# manual

---


## Body scale

Item: 48834-48835

EN

NO

SV

[www.rubicson.com](http://www.rubicson.com)  
Box 50435, Malmö, Sweden  
2018-01-16

**Rubicson®**

## Specifications

Batteries:	3x AAA
Capacity:	180 kg
Size:	303x24x303 mm

## Start

Insert 3x AAA batteries in the battery hatch. Mind the polarity.

## Standard settings

- Place the scale on a flat surface to ensure accuracy.
- Step on the scale with one foot, it will turn on and display 0.0, remove your foot from the scale. Wait for the scale to turn off. This process needs to be repeated if the scale is moved.
- The scale is now ready for use, step on the scale and stand still while your weight is measured.
- The scale automatically turns off if no weighing is done after 10 seconds.
- You can change the weight unit by pressing the plastic tab labeled UNIT on the bottom of the scale. Switch between kilograms (kg), pounds (lb) and stones (st).

## Troubleshooting

- When the display says "LO", insert new batteries.
- When the display says "Err", it means that the maximum weight has been exceeded.

## Warning

Make sure that the scale is placed on a flat surface to ensure accuracy.

Always weigh yourself without clothing or footwear, and before meals. The first thing in the morning is the best time.

Keep the scale in a dry place to protect the electronic components.

Do not store the scale in an upright position as this will drain the batteries.

Avoid overloading the scale; otherwise, permanent damage may occur.

Remove the batteries if the scales are not going to be used for a long period of time.

Treat your scale with care as it is a precision instrument, do not drop it or jump up or down on it

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